

Vathin Single-Use Ureteral Access Sheath

Unlocking Pathways to
Successful Procedures

Hunan Vathin Medical Instrument Co.,Ltd.


1/F, Building 12, Innovation and Entrepreneurship Service Center No. 9
Chuanqi West Road, Jiuhe Economic Development Zone,
411100 Xiangtan, Hunan, China

 www.vathin.com  marketing@vathin.com

MD24020

Vathin Single-Use Ureteral Access Sheath

The Ultimate Solution for Smooth Scope Manipulation and Repeated Passage



| Model | Sheath Specification | Sheath Effective Length(mm) | Sheath Outer Diameter (mm) | Sheath Inner Diameter (mm) | Dilator Length (mm) | Dilator Outer Diameter (mm) | Dilator Inner Diameter (mm) | Tip Bending Range |
|----------|----------------------|-----------------------------|----------------------------|----------------------------|---------------------|-----------------------------|-----------------------------|-------------------|
| UA-S1800 | 10/12Fr | 350 | 4.13 | 3.33 | 448 | 3.3 | 1.05 | >200° |
| UA-S1801 | 10/12Fr | 450 | 4.13 | 3.33 | 548 | 3.3 | | |
| UA-S1A00 | 11/13Fr | 350 | 4.4 | 3.7 | 448 | 3.65 | | |
| UA-S1A01 | 11/13Fr | 450 | 4.4 | 3.7 | 548 | 3.65 | | |
| UA-S1B00 | 12/14Fr | 350 | 4.83 | 3.97 | 448 | 3.95 | | |
| UA-S1B01 | 12/14Fr | 450 | 4.83 | 3.97 | 548 | 3.95 | | |
| UA-S1C00 | 13/15Fr | 350 | 5.1 | 4.3 | 448 | 4.25 | | |
| UA-S1C01 | 13/15Fr | 450 | 5.1 | 4.3 | 548 | 4.25 | | |
| UA-S1D00 | 14/16Fr | 350 | 5.53 | 4.64 | 448 | 4.6 | | |
| UA-S1D01 | 14/16Fr | 450 | 5.53 | 4.64 | 548 | 4.6 | | |