

Vathin Single-use Ureteral Access Sheath

Unlocking Pathways to
Successful Procedures

Hunan Vathin Medical Instrument Co.,Ltd.


1/F, Building 12, Innovation and Entrepreneurship Service Center No. 9
Chuanqi West Road, Jiuhe Economic Development Zone,
411100 Xiangtan, Hunan, China

 www.vathin.com

 marketing@vathin.com

Negative Pressure For Clearer Vision & Faster Stone Removal

Enhanced Clarity and Efficiency



Model	Sheath Specification	Sheath Effective Length(mm)	Sheath Outer Diameter (mm)	Sheath Inner Diameter (mm)	Dilator Length (mm)	Dilator Outer Diameter (mm)	Dilator Inner Diameter (mm)	Tip Bending Range
UA-S1800	10/12Fr	350	4.13	3.33	448	3.3	1.05	>200°
UA-S1801	10/12Fr	450	4.13	3.33	548	3.3		
UA-S1A00	11/13Fr	350	4.4	3.7	448	3.65		
UA-S1A01	11/13Fr	450	4.4	3.7	548	3.65		
UA-S1B00	12/14Fr	350	4.83	3.97	448	3.95		
UA-S1B01	12/14Fr	450	4.83	3.97	548	3.95		